

# ZNF586 Antibody



PACO31400

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF586

**Uniprot**

Q9NXT0

**Synonyms:**

Zinc finger protein 586, ZNF586

**Immunogen:**


Recombinant Human Zinc finger protein 586 protein (1-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO31400 at dilution of 1:100.