

ZNF610 Antibody



PACO31416

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF610

Uniprot

Q8N9Z0

Synonyms:

Zinc finger protein 610, ZNF610

Immunogen:

Recombinant Human Zinc finger protein 610 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF610 antibody at 7 μ g/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 54, 49 kDa

Observed band size: 54 kDa