

ZNF614 Antibody



PACO31420

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF614

Uniprot

Q8N883

Synonyms:

Zinc finger protein 614, ZNF614

Immunogen:

Recombinant Human Zinc finger protein 614 protein (1-198AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF614 antibody at 6 μ g/ml

Lane 1: Jurkat whole cell lysate

Lane 2: HepG2 whole cell lysate


Lane 3: A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 68, 23 kDa

Observed band size: 68 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31420 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO31420 at dilution of 1:100.