

ZNF765 Antibody



PACO31488

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF765

Uniprot

Q7L2R6

Synonyms:

Zinc finger protein 765, ZNF765

Immunogen:

Recombinant Human Zinc finger protein 765 protein (1-60AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A