

ZNF77 Antibody



PACO31492

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF77

Uniprot

Q15935

Synonyms:

Zinc finger protein 77 (ZNFpT1), ZNF77


Immunogen:

Recombinant Human Zinc finger protein 77 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31492 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO31492 at dilution of 1:100.