

ZNF791 Antibody



PACO31500

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF791

Uniprot

Q3KP31

Synonyms:

Zinc finger protein 791, ZNF791

Immunogen:

Recombinant Human Zinc finger protein 791 protein (1-250AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO31500 at dilution of 1:100.



Immunohistochemistry analysis of human colon cancer using PACO31500 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO31500 at dilution of 1:100.