

ZNF821 Antibody



PACO31512

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF821

Uniprot

O75541

Synonyms:

Zinc finger protein 821, ZNF821

Immunogen:

Recombinant Human Zinc finger protein 821 protein (1-370AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF821 antibody at 8 μ g/ml + MCF-7 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 47, 43 kDa

Observed band size: 47 kDa

.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO31512 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO31512 at dilution of 1:100.