

# ZNF98 Antibody



PACO31524

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:50-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF98

**Uniprot**

A6NK75

**Synonyms:**

Zinc finger protein 98 (Zinc finger protein 739) (Zinc finger protein F7175), ZNF98, ZNF739

**Immunogen:**


Recombinant Human Zinc finger protein 98 protein (1-135AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunofluorescent analysis of HeLa cells using PACO31524 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).