

PACO31556

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Lysozymes have primarily a bacteriolytic function; those in tissues and body fluids are associated with the monocyte-macrophage system and enhance the activity of immunoagents.

Gene ID:

LYZ

Uniprot

P61626

Synonyms:

Lysozyme C (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase C), LYZ, LZM


Immunogen:

Recombinant Human Lysozyme C protein (19-148AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO31556 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO31556 at dilution of 1:100.