

# STARD7 Antibody



PACO31672

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

May play a protective role in mucosal tissues by preventing exaggerated allergic responses.

**Gene ID:**

STARD7

**Uniprot**

Q9NQZ5

**Synonyms:**

StAR-related lipid transfer protein 7, mitochondrial (Gestational trophoblastic tumor protein 1) (START domain-containing protein 7) (StARD7), STARD7, GTT1

**Immunogen:**


Recombinant Human StAR-related lipid transfer protein 7, mitochondrial protein (61-307AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western blot

All lanes: STARD7 antibody at 2 $\mu$ g/ml + PC-3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa