

PACO31884

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. Down-regulates activation of NF-kappa-B.

Gene ID:

COMMD4

Uniprot

Q9H0A8

Synonyms:

COMM domain-containing protein 4, COMMD4


Immunogen:

Recombinant Human COMM domain-containing protein 4 protein (1-195AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO31884 at dilution of 1:100.