

CTDNEP1 Antibody



PACO31964

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CTDNEP1

Uniprot

O95476

Synonyms:

CTD nuclear envelope phosphatase 1 (EC 3.1.3.16) (Serine/threonine-protein phosphatase dullard), CTDNEP1, DULLARD


Immunogen:

Recombinant Human CTD nuclear envelope phosphatase 1 protein (3-244AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31964 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO31964 at dilution of 1:100.