

PACO32132

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

EGF stimulates the growth of various epidermal and epithelial tissues in vivo and in vitro and of some fibroblasts in cell culture. Magnesiotropic hormone that stimulates magnesium reabsorption in the renal distal convoluted tubule via engagement of EGFR and activation of the magnesium channel TRPM6.

Gene ID:

EGF

Uniprot

P01133

Synonyms:

Pro-epidermal growth factor (EGF) [Cleaved into: Epidermal growth factor (Urogastrone)], EGF

Immunogen:

Recombinant Human Pro-epidermal growth factor protein (977-1023AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A