

PACO32252

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

LCN2

Uniprot

P80188

Synonyms:

Neutrophil gelatinase-associated lipocalin (NGAL) (25 kDa alpha-2-microglobulin-related subunit of MMP-9) (Lipocalin-2) (Oncogene 24p3) (Siderocalin LCN2) (p25), LCN2, HNL NGAL


Immunogen:

Recombinant Human NGAL protein (1-198AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: LCN2 antibody at 2 μ g/ml

Lane 1: Human positive serum at 1: 10

Lane 2: Human positive serum at 1: 100

Lane 3: Human positive serum at 1: 500

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa