

IL1RL1 Antibody



PACO32348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for interleukin-33 (IL-33), its stimulation recruits MYD88, IRAK1, IRAK4, and TRAF6, followed by phosphorylation of MAPK3/ERK1 and/or MAPK1/ERK2, MAPK14, and MAPK8. Possibly involved in helper T-cell function.

Gene ID:

IL1RL1

Uniprot

Q01638

Synonyms:

Interleukin-1 receptor-like 1 (Protein ST2), IL1RL1, DER4 ST2 T1

Immunogen:

Recombinant Human Interleukin-1 receptor-like 1 protein (19-328AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A