

TGFB2 Antibody, FITC conjugated



PACO32374

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

TGF-beta 2 has suppressive effects on interleukin-2 dependent T-cell growth.

Gene ID:

TGFB2

Uniprot

P61812

Synonyms:

Transforming growth factor beta-2 (TGF-beta-2) (BSC-1 cell growth inhibitor) (Cetermin) (Glioblastoma-derived T-cell suppressor factor) (G-TSF) (Polyergin) [Cleaved into: Latency-associated peptide (LAP)], TGFB2

Immunogen:

Recombinant Human Transforming growth factor beta-2 proprotein protein (304-413AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A