

THBS1 Antibody, FITC conjugated



PACO32378

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

THBS1

Uniprot

P07996

Synonyms:

Thrombospondin-1, THBS1, TSP TSP1

Immunogen:

Recombinant Human Thrombospondin-1 protein (19-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A