

TIMP4 Antibody, Biotin conjugated



PACO32387

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by binding to their catalytic zinc cofactor. Known to act on MMP-1, MMP-2, MMP-3, MMP-7 and MMP-9.

Gene ID:

TIMP4

Uniprot

Q99727

Synonyms:

Metalloproteinase inhibitor 4 (Tissue inhibitor of metalloproteinases 4) (TIMP-4), TIMP4

Immunogen:

Recombinant Human Metalloproteinase inhibitor 4 protein (31-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A