

L1CAM Antibody, Biotin conjugated



PACO32395

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cell adhesion molecule with an important role in the development of the nervous system. Involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc. Binds to axonin on neurons.

Gene ID:

L1CAM

Uniprot

P32004

Synonyms:

Neural cell adhesion molecule L1 (N-CAM-L1) (NCAM-L1) (CD antigen CD171), L1CAM, CAML1 MIC5

Immunogen:

Recombinant Human Neural cell adhesion molecule L1 protein (995-1119AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A