

blaNDM-1 Antibody, Biotin conjugated



PACO32435

Product Information

Size:

50ug

Reactivity:

Klebsiella pneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Confers resistance to many beta-lactam antibiotics, including some carbapenems. Does not confer resistance to the polymixin colistin or the fluoroquinolone ciprofloxacin.

Gene ID:

blaNDM-1

Uniprot

C7C422

Synonyms:

Metallo-beta-lactamase type 2 (EC 3.5.2.6) (B2 metallo-beta-lactamase) (Beta-lactamase type II) (Metallo-beta-lactamase NDM-1) (Metallo-beta-lactamase type II) (New Delhi metallo-beta-lactamase-1) (NDM-1), blaNDM-1

Immunogen:

Recombinant Klebsiella pneumoniae Metallo-beta-lactamase type 2 protein (29-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A