

Thpo Antibody



PACO32440

Product Information

Size:

50ug

Reactivity:

Rat, Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Lineage-specific cytokine affecting the proliferation and maturation of megakaryocytes from their committed progenitor cells. It acts at a late stage of megakaryocyte development. It may be the major physiological regulator of circulating platelets.

Gene ID:

Thpo

Uniprot

P49745

Synonyms:

Thrombopoietin, Thpo

Immunogen:

Recombinant Rat Thrombopoietin protein (22-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Thpo antibody at 4 μ g/ml

Lane 1: Mouse liver tissue

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

.