

# FABP3 Antibody, FITC conjugated



PACO32458

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FABP3

**Uniprot**

P05413

**Synonyms:**

Fatty acid, binding protein, heart (Fatty acid, binding protein 3) (Heart-type fatty acid, binding protein) (H-FABP) (Mammary-derived growth inhibitor) (MDGI) (Muscle fatty acid, binding protein) (M-FABP), FABP3, FABP11 MDGI

**Immunogen:**

Recombinant Human Fatty acid, binding protein, heart protein (2-133AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A