

PACO32512

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Gene ID:

TNNI3

Uniprot

P19429

Synonyms:

Troponin I, cardiac muscle (Cardiac troponin I), TNNI3, TNNC1

Immunogen:

Recombinant Human Troponin I, cardiac muscle protein (2-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A