

CALM1 Antibody, HRP conjugated



PACO32581

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

CALM2 is the second calmodulin gene. Although the CALM1, CALM2, and CALM3 calmodulin proteins are identical, at the nucleotide level they share only about 80% identity within their coding regions, and they contain no significant homology within their noncoding regions.

Gene ID:

CALM1

Uniprot

P62158

Synonyms:

Calmodulin-1, CALM1, CALM CAM CAM1

Immunogen:

Recombinant Human Calmodulin protein (2-149AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A