

IL18 Antibody, HRP conjugated



PACO32597

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IL18

Uniprot

Q14116

Synonyms:

Interleukin-18 (IL-18) (Iboctadecin) (Interferon gamma-inducing factor) (IFN-gamma-inducing factor) (Interleukin-1 gamma) (IL-1 gamma), IL18, IGIF IL1F4

Immunogen:

Recombinant Human Interleukin-18 protein (37-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A