

# IL18 Antibody, FITC conjugated



PACO32598

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IL18

**Uniprot**

Q14116

**Synonyms:**

Interleukin-18 (IL-18) (Iboctadecin) (Interferon gamma-inducing factor) (IFN-gamma-inducing factor) (Interleukin-1 gamma) (IL-1 gamma), IL18, IGIF IL1F4

**Immunogen:**

Recombinant Human Interleukin-18 protein (37-193AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A