

IGF1 Antibody, HRP conjugated



PACO32605

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IGF1

Uniprot

P05019

Synonyms:

Insulin-like growth factor I (IGF-I) (Mechano growth factor) (MGF) (Somatomedin-C), IGF1, IBP1

Immunogen:

Recombinant Human Insulin-like growth factor I protein (49-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A