

Ccl4 Antibody, HRP conjugated



PACO32657

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Monokine with inflammatory and chemokinetic properties.

Gene ID:

Ccl4

Uniprot

P14097

Synonyms:

C-C motif chemokine 4 (Immune activation protein 2) (ACT-2) (ACT2) (Macrophage inflammatory protein 1-beta) (MIP-1-beta) (Protein H400) (SIS-gamma) (Small-inducible cytokine A4), Ccl4, Mip1b Scya4

Immunogen:

Recombinant Mouse C-C motif chemokine 4 protein (24-92AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A