

# Cxcl1 Antibody, HRP conjugated



PACO32673

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Has chemotactic activity for neutrophils. Contributes to neutrophil activation during inflammation. By similarity. Hematopoietic chemokine, which, in vitro, suppresses hematopoietic progenitor cell proliferation. KC(5-72) shows a highly enhanced hematopoietic activity.

**Gene ID:**

Cxcl1

**Uniprot**

P12850

**Synonyms:**

Growth-regulated alpha protein (C-X-C motif chemokine 1) (Platelet-derived growth factor-inducible protein KC) (Secretory protein N51) [Cleaved into: KC(5-72) (Hematopoietic synergistic factor) (HSF) (KC-T)], Cxcl1, Gro Gro1 Mgsa Scyb1

**Immunogen:**

Recombinant Mouse Growth-regulated &alpha; protein (29-96AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A