

Cxcl1 Antibody, Biotin conjugated



PACO32675

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has chemotactic activity for neutrophils. Contributes to neutrophil activation during inflammation. By similarity. Hematopoietic chemokine, which, in vitro, suppresses hematopoietic progenitor cell proliferation. KC(5-72) shows a highly enhanced hematopoietic activity.

Gene ID:

Cxcl1

Uniprot

P12850

Synonyms:

Growth-regulated alpha protein (C-X-C motif chemokine 1) (Platelet-derived growth factor-inducible protein KC) (Secretory protein N51) [Cleaved into: KC(5-72) (Hematopoietic synergistic factor) (HSF) (KC-T)], Cxcl1, Gro Gro1 Mgsa Scyb1

Immunogen:

Recombinant Mouse Growth-regulated α protein (29-96AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A