

Cxcl11 Antibody, Biotin conjugated



PACO32691

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Cxcl11

Uniprot

Q9JHH5

Synonyms:

C-X-C motif chemokine 11 (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11), Cxcl11, Scyb11

Immunogen:

Recombinant Mouse C-X-C motif chemokine 11 protein (22-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A