

Cxcl14 Antibody, Biotin conjugated



PACO32695

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Chemotactic for CESS B-cells and THP-1 monocytes, but not T-cells.

Gene ID:

Cxcl14

Uniprot

Q9WUQ5

Synonyms:

C-X-C motif chemokine 14 (B-cell and monocyte-activating chemokine) (Chemokine BRAK) (Kidney-expressed chemokine CXC) (MIP-2G) (Small-inducible cytokine B14), Cxcl14, Bmac Kec Ks1 Mip2g Scyb14

Immunogen:

Recombinant Mouse C-X-C motif chemokine 14 protein (23-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A