

# LGALS7 Antibody, Biotin conjugated



PACO32731

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Could be involved in cell-cell and/or cell-matrix interactions necessary for normal growth control. Pro-apoptotic protein that functions intracellularly upstream of JNK activation and cytochrome c release.

**Gene ID:**

LGALS7

**Uniprot**

P47929

**Synonyms:**

Galectin-7 (Gal-7) (HKL-14) (PI7) (p53-induced gene 1 protein), LGALS7; LGALS7B, PIG1;

**Immunogen:**

Recombinant Human Galectin-7 protein (2-136AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A