

# PTH Antibody, HRP conjugated



PACO32753

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-14C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

**Gene ID:**

PTH

**Uniprot**

P01270

**Synonyms:**

Parathyroid hormone (PTH) (Parathormone) (Parathyrin), PTH

**Immunogen:**

Recombinant Human Parathyroid hormone protein (32-115AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A