

PTH Antibody, Biotin conjugated



PACO32755

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-14C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

Gene ID:

PTH

Uniprot

P01270

Synonyms:

Parathyroid hormone (PTH) (Parathormone) (Parathyrin), PTH

Immunogen:

Recombinant Human Parathyroid hormone protein (32-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A