

Lgals7 Antibody, Biotin conjugated



PACO32775

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Could be involved in cell-cell and/or cell-matrix interactions necessary for normal growth control. Pro-apoptotic protein that functions intracellularly upstream of JNK activation and cytochrome c release.

Gene ID:

Lgals7

Uniprot

O54974

Synonyms:

Galectin-7 (Gal-7), Lgals7

Immunogen:

Recombinant Mouse Galectin-7 protein (2-136AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A