

STC1 Antibody, Biotin conjugated



PACO32779

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

Gene ID:

STC1

Uniprot

P52823

Synonyms:

Stanniocalcin-1 (STC-1), STC1, STC

Immunogen:

Recombinant Human Stanniocalcin-1 protein (39-247AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A