

NTF4 Antibody, Biotin conjugated



PACO32787

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Target-derived survival factor for peripheral sensory sympathetic neurons.

Gene ID:

NTF4

Uniprot

P34130

Synonyms:

Neurotrophin-4 (NT-4) (Neurotrophin-5) (NT-5) (Neurotrophic factor 4), NTF4, NTF5

Immunogen:

Recombinant Human Neurotrophin-4 protein (82-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A