

AFP Antibody, Biotin conjugated



PACO32791

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds copper, nickel, and fatty acid, as well as, and bilirubin less well than, serum albumin. Only a small percentage (less than 2%) of the human AFP shows estrogen-binding properties.

Gene ID:

AFP

Uniprot

P02771

Synonyms:

Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), AFP, HPAFP

Immunogen:

Recombinant Human Alpha-fetoprotein protein (19-609AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A