

TNNT2 Antibody, HRP conjugated



PACO32809

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Gene ID:

TNNT2

Uniprot

P45379

Synonyms:

Troponin T, cardiac muscle (TnTc) (Cardiac muscle troponin T) (cTnT), TNNT2

Immunogen:

Recombinant Human Troponin T, cardiac muscle protein (138-298AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A