

DCN Antibody, HRP conjugated



PACO32833

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May affect the rate of fibrils formation.

Gene ID:

DCN

Uniprot

P07585

Synonyms:

Decorin (Bone proteoglycan II) (PG-S2) (PG40), DCN, SLRR1B

Immunogen:

Recombinant Human Decorin protein (20-359AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A