

# DCN Antibody, FITC conjugated



PACO32834

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May affect the rate of fibrils formation.

**Gene ID:**

DCN

**Uniprot**

P07585

**Synonyms:**

Decorin (Bone proteoglycan II) (PG-S2) (PG40), DCN, SLRR1B

**Immunogen:**

Recombinant Human Decorin protein (20-359AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A