

GSTP1 Antibody, HRP conjugated



PACO32837

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration.

Gene ID:

GSTP1

Uniprot

P09211

Synonyms:

Glutathione S-transferase P (EC 2.5.1.18) (GST class-pi) (GSTP1-1), GSTP1, FAEES3 GST3

Immunogen:

Recombinant Human Glutathione S-transferase P protein (2-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A