

# GSTP1 Antibody, FITC conjugated



PACO32838

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration.

**Gene ID:**

GSTP1

**Uniprot**

P09211

**Synonyms:**

Glutathione S-transferase P (EC 2.5.1.18) (GST class-pi) (GSTP1-1), GSTP1, FAEES3 GST3

**Immunogen:**

Recombinant Human Glutathione S-transferase P protein (2-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A