

GKN1 Antibody, Biotin conjugated



PACO32867

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has mitogenic activity and may be involved in maintaining the integrity of the gastric mucosal epithelium.

Gene ID:

GKN1

Uniprot

Q9NS71

Synonyms:

Gastroke-1 (18 kDa antrum mucosa protein) (AMP-18) (Protein CA11), GKN1, AMP18 CA11

Immunogen:

Recombinant Human Gastroke-1 protein (49-199AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A