

BMP6 Antibody, FITC conjugated



PACO32894

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Induces cartilage and bone formation.

Gene ID:

BMP6

Uniprot

P22004

Synonyms:

Bone morphogenetic protein 6 (BMP-6) (VG-1-related protein) (VG-1-R) (VGR-1), BMP6, VGR

Immunogen:

Recombinant Human Bone morphogenetic protein 6 protein (382-513AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A