

CDH3 Antibody, HRP conjugated



PACO32897

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Gene ID:

CDH3

Uniprot

P22223

Synonyms:

Cadherin-3 (Placental cadherin) (P-cadherin), CDH3, CDHP

Immunogen:

Recombinant Human Cadherin-3 protein (117-360AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A