

RARRES2 Antibody, FITC conjugated



PACO32910

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RARRES2

Uniprot

Q99969

Synonyms:

Retinoic acid, receptor responder protein 2 (Chemerin) (RAR-responsive protein TIG2) (Tazarotene-induced gene 2 protein), RARRES2, TIG2

Immunogen:

Recombinant Human Retinoic acid, receptor responder protein 2 protein (21-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A