

HSPB7 Antibody



PACO32936

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

HSPB7

Uniprot

Q9UBY9

Synonyms:

Heat shock protein beta-7 (HspB7) (Cardiovascular heat shock protein) (cvHsp), HSPB7, CVHSP


Immunogen:

Recombinant Human Heat shock protein beta-7 protein (1-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO32936 at dilution of 1:100.