

HSPA12B Antibody



PACO32940

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

HSPA12B

Uniprot

Q96MM6

Synonyms:

Heat shock 70 kDa protein 12B, HSPA12B, C20orf60


Immunogen:

Recombinant Human Heat shock 70 kDa protein 12B protein (201-500AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO32940 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO32940 at dilution of 1:100.